

**SHENTEK**

# **ResHCD™ Human DNA Quantitation Kit User Guide**

Version: A/0

For Research Use Only

Product No.: 1101119

Reagents for 100 Reactions

**Huzhou Shenke Biotechnology Co., Ltd**

*(IMPORTANT: Please read this document carefully before experiment.)*

## 1. Product information

### ■ Product description

The ResHCD™ Human DNA Quantitation Kit is used for quantitation of residual Human host cell DNA in a variety of biopharmaceutical products. By combining quantitative PCR, the kit ensures a rapid and specific quantitation of residual DNA at the femtogram level. For extraction information, please refer to the SHENTEK® Residual Host Cell DNA Sample Preparation Kit (Fast) User Guide (Product No. 1104201).

The kit is for Research Use Only and is not for use in diagnostic procedures or for therapeutic purposes in humans or animals.

### ■ Kit components and storage

*WARNING: Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, mask, clothing and gloves.*

Table 1. Kit components and storage

| Reagent                   | Part No. | Quantity         | Storage                         |
|---------------------------|----------|------------------|---------------------------------|
| qPCR Reaction MIX-2       | NNB029   | 850 µL × 2 tubes | -20°C,<br>protect from<br>light |
| Human Primer&Probe MIX    | NNC007   | 300 µL × 1 tube  |                                 |
| IPC MIX-4                 | NNC143   | 200 µL × 1 tube  |                                 |
| DNA Dilution Buffer (DDB) | NND001   | 1.5 mL × 3 tubes | -20°C                           |
| Human DNA Control         | NNA003   | 50 µL × 1 tube   |                                 |

*Note: According to USP <659>, items with recommended storage temperatures not exceeding –20°C shall be stored within ± 10°C (–30°C to –10°C). Short-term storage below –30°C is acceptable when supported by stability data.*

The kit components can be stored under the appropriate conditions for up to 24 months. Please check the expiration date on the labels.

## ■ Applied instruments, including but not limited to the following

- SHENTEK-96H Real-Time PCR System
- ABI 7500 Real-Time PCR System
- LightCycler 480 II Real-Time PCR System


## ■ Required materials not included in the kit

- Low retention, RNase/DNase-free, sterile microcentrifuge tubes
- 96-well qPCR plates with sealing film, or 8-tube PCR strips with caps
- Low retention filter tips: 1000  $\mu$ L, 100  $\mu$ L and 10  $\mu$ L

## ■ Related equipments

- Real-Time PCR System
- Vortex mixer
- Benchtop microcentrifuge
- Pipettes: 1000  $\mu$ L, 100  $\mu$ L and 10  $\mu$ L
- Microplate and microtube shaker

## ■ Workflow



## 2. Methods

### ■ Experiment preparation

1. Wear appropriate protective eyewear, mask, clothing and gloves.
2. Irradiate the tabletop, pipettes and tubes with UV for 30 minutes, and disinfect with 75% ethanol.
3. Thaw the kit completely at 2-8°C or melt on ice, vortex and spin briefly.

### ■ DNA Control serial dilutions for the standard curve

Please check the concentration labeled on the tube containing the Human DNA Control prior to dilution.

Prepare a series of Human DNA Control solution with DNA Dilution Buffer (DDB) as follows:

1. Thaw Human DNA Control and DNA Dilution Buffer completely at 2-8°C or melt on ice. Vortex to mix well and centrifuge briefly, repeat 3 times.
2. Label six 1.5 mL microcentrifuge tubes as ST0, ST1, ST2, ST3, ST4 and ST5 respectively.
3. Dilute the Human DNA Control to 3000 pg/μL with DDB in the ST0 tube.

Calculate the volume of DDB to prepare the ST0:

$$\frac{\text{DNA Control conc. (A)} \times 1000 \text{ pg/ng} \times \text{Volume of DNA Control (B)}}{3000 \text{ pg/}\mu\text{L}} - \text{Volume of DNA Control (B)}$$

For example:

The concentration on the label of the DNA Control is 30.9 ng/μL (A), pipette 10 μL (B) of the DNA Control to the ST0 tube. Add DDB in the volume shown in the calculation below to reach a final concentration of 3000 pg/μL.

$$\frac{30.9 \text{ ng/}\mu\text{L} \times 1000 \text{ pg/ng} \times 10 \mu\text{L}}{3000 \text{ pg/}\mu\text{L}} - 10 \mu\text{L} = 93 \mu\text{L}$$

4. Vortex to mix well and centrifuge briefly the ST0 tube for 3-5 seconds in microcentrifuge, and repeat 3 times to mix thoroughly.
5. Add 90 μL DDB to each tube of ST1, ST2, ST3, ST4 and ST5.
6. Perform the serial dilutions according to Table 2:

Table 2. Dilution for Human DNA Control

| Serial dilution tube | Dilution                        | Conc. (pg/ $\mu$ L) |
|----------------------|---------------------------------|---------------------|
| ST0                  | Dilute the DNA control with DDB | 3000                |
| ST1                  | 10 $\mu$ L ST0 + 90 $\mu$ L DDB | 300                 |
| ST2                  | 10 $\mu$ L ST1 + 90 $\mu$ L DDB | 30                  |
| ST3                  | 10 $\mu$ L ST2 + 90 $\mu$ L DDB | 3                   |
| ST4                  | 10 $\mu$ L ST3 + 90 $\mu$ L DDB | 0.3                 |
| ST5                  | 10 $\mu$ L ST4 + 90 $\mu$ L DDB | 0.03                |

*Note: The remaining unused DDB need to be stored at 2-8°C. If the solution is cloudy or contains precipitates, heat at 37°C until it clear.*

*Note: At least five concentration values should be included when creating a standard curve. To select appropriate sample dilutions, we recommend performing method validation before sample testing.*

## ■ Sample preparation

### ➤ Preparation for Extraction Recovery Control (ERC) samples

According to the Human DNA spike concentration in ERC samples (Take the samples containing 30 pg of Human DNA as an example), the specific preparation procedure is as follows:

1. Aliquot 100  $\mu$ L of the test sample to a new 1.5 mL microcentrifuge tube.
2. Add 10  $\mu$ L of ST3 solution and mix thoroughly, label it as “ERC sample.

### ➤ Preparation for Negative Control Samples (NCS)

Add 100  $\mu$ L of DDB to a new 1.5 mL microcentrifuge tube, and label it as NCS.

The ERC sample and NCS should be processed in the same procedures as test sample preparation before testing.

## ■ qPCR Reaction MIX preparation

1. Determine the number of reaction wells based on your selected standard curve, with the number of test samples and control samples. Generally, triplicates are tested for each sample.

Number of reaction wells = (standard curve of 5 concentration gradients + 1 NTC + 1 NCS + test samples )  $\times$  3.

2. Prepare qPCR Reaction MIX according to Table 3 and Table 4

3. Refer to 96-well plate layouts shown in Table 5.

Table 3. qPCR MIX Preparation

| Reagents               | Volume/reaction | Volume for 30 reactions<br>(includes 10% overage) |
|------------------------|-----------------|---------------------------------------------------|
| qPCR Reaction MIX-2    | 15 $\mu$ L      | 495 $\mu$ L                                       |
| Human Primer&Probe MIX | 3 $\mu$ L       | 99 $\mu$ L                                        |
| IPC MIX-4              | 2 $\mu$ L       | 66 $\mu$ L                                        |
| Total volume           | 20 $\mu$ L      | 660 $\mu$ L                                       |

Table 4. qPCR Reaction MIX Preparation

|                |                                                       |
|----------------|-------------------------------------------------------|
| Standard curve | 20 $\mu$ L qPCR MIX + 10 $\mu$ L ST1/ST2/ST3/ST4/ST5  |
| NTC            | 20 $\mu$ L qPCR MIX + 10 $\mu$ L DDB                  |
| NCS            | 20 $\mu$ L qPCR MIX + 10 $\mu$ L purified NCS         |
| Test sample    | 20 $\mu$ L qPCR MIX + 10 $\mu$ L purified test sample |
| ERC sample     | 20 $\mu$ L qPCR MIX + 10 $\mu$ L purified ERC sample  |

Table 5. Example of 96-well plate layouts

|     |   |    |    |    |           |           |           |   |     |     |     |   |
|-----|---|----|----|----|-----------|-----------|-----------|---|-----|-----|-----|---|
| NTC |   | S1 | S1 | S1 | S1<br>ERC | S1<br>ERC | S1<br>ERC |   | ST5 | ST5 | ST5 | A |
| NTC |   | S2 | S2 | S2 | S2<br>ERC | S2<br>ERC | S2<br>ERC |   | ST4 | ST4 | ST4 | B |
| NTC |   | S3 | S3 | S3 | S3<br>ERC | S3<br>ERC | S3<br>ERC |   | ST3 | ST3 | ST3 | C |
|     |   | S4 | S4 | S4 | S4<br>ERC | S4<br>ERC | S4<br>ERC |   | ST2 | ST2 | ST2 | D |
| NCS |   | S5 | S5 | S5 | S5<br>ERC | S5<br>ERC | S5<br>ERC |   | ST1 | ST1 | ST1 | E |
| NCS |   |    |    |    |           |           |           |   |     |     |     | F |
| NCS |   |    |    |    |           |           |           |   |     |     |     | G |
|     |   |    |    |    |           |           |           |   |     |     |     | H |
| 1   | 2 | 3  | 4  | 5  | 6         | 7         | 8         | 9 | 10  | 11  | 12  |   |

*Note: This example represents the assay for a standard curve with 5 concentration gradients (ST1 to ST5), 1 NTC, 1 NCS, and 5 test samples (S1 to S5), and 5 ERC sample (S1 ERC to S5 ERC), and 3 replicates for each sample.*

*Note: In specific testing, the plate layouts for sample loading can be adjusted based on the sample quantity.*

4. Seal the 96-well plate with sealing film. Mix contents of the wells in microplate shaker, then spin down the plate for 10 seconds in microcentrifuge and place it in the qPCR instrument.

## ■ qPCR program setting

*Note: The following instructions apply only to the ABI7500 Real-Time PCR System with SDS v1.5.1. If you use a different instrument or software, refer to the applicable instrument or software documentation.*

1. Create a new document, then in the Assay drop-down list, select **Standard Curve (Absolute Quantitation)**.
2. In the Run Mode drop-down list, select **Standard 7500**, then click **Next**.
3. Click **New Detector**:
  - a. Enter **Human-DNA** in the Name field.
  - b. Select **FAM** in the Reporter Dye drop-down list and select **(none)** in the

- Quencher Dye drop-down list, then click **OK**.
- c. Select a color for the detector, then click **Create Another**.
4. Click **New Detector**:
- a. Enter **IPC** in the Name field.
  - b. Select **VIC** in the Reporter Dye drop-down list and select **(none)** in the Quencher Dye drop-down list, then click **OK**.
  - c. Select a color for the detector, then click **OK**.
  - d. Select the detectors, then click **Add** to add the detectors to the document.
5. Select **ROX** as the passive reference dye, then Click **Next**.
6. Select the applicable set of wells for the samples, then select Human-DNA and IPC detector for each well.
7. Select **Finish**, and then set thermal-cycling conditions:
- a. Choose the thermal cycling reaction volume to 30  $\mu$ L.
  - b. Set up the program as following:

Table 6. qPCR running temperature and time

| Step                | Temp. | Time(mm:sec) | Cycles |
|---------------------|-------|--------------|--------|
| Activation          | 95°C  | 10:00        | 1      |
| Denaturation        | 95°C  | 00:15        | 40     |
| Annealing/extension | 60°C* | 1:00         |        |

\*Instrument will read the fluorescence signal during this step.

8. Save the document, then click **Start** to start the real-time qPCR run.

## ■ Results analysis

1. Select **Set up** tab, then set tasks for each sample type by clicking on the Task Column drop-down list:
  - a. NTC: target DNA detector task = **NTC**
  - b. NCS, test samples, and ERC wells: target DNA detector task = **Unknown**
2. Set up the standard curve as shown in the following table:

Table 7. Settings for Standard curve

| Tube label | Task     | Conc. (pg/ $\mu$ L) |
|------------|----------|---------------------|
| ST1        | Standard | 300                 |
| ST2        | Standard | 30                  |
| ST3        | Standard | 3                   |
| ST4        | Standard | 0.3                 |
| ST5        | Standard | 0.03                |

3. Select the **Results** tab, then select Amplification Plot.
4. In the Data drop-down list, select **Delta Rn vs Cycle**.
5. In the Analysis Settings window, enter the following settings:
  - a. Select **Manual Ct**.
  - b. In the Threshold field, enter value 0.05 of FAM; value 0.1 of VIC.
  - c. Select **Automatic Baseline**.
6. Click the button  in the toolbar, then wait the plate analyzing.
7. Select the **Result** tab> >**Standard curve** tab, then verify the Slope, Intercept and  $R^2$  values.
8. Select the Report tab, then achieve the mean quantity and standard deviation for each sample.
9. Select **File > > Export > > Results**. In the Save as type drop-down list, select **Results Export Files**, then click **Save**.
10. The recovery of ERC samples is calculated based on the value of the test samples and the ERC samples. The recovery should be between 50% and 150%.
11. The Ct value of NCS should be higher than the mean Ct value of the lowest concentration in the standard curve. If the validated quantitation limit (QL) concentration is lower than the lowest concentration in the standard curve, the value of the NCS should be below the QL.
12. The detection value of NTC should be no more than 6.00 fg/ $\mu$ L, or set specific criteria by your own method validation, meanwhile shows normal

amplification curve in the VIC signal channel.

*Note: The parameter settings listed above are generally recommended for result analysis and could be adjusted depending on the specific instrument model and software version. In most cases, the instrument performs the interpretation automatically.*

Effective date: 20 Apr. 2026

## Support & Contact

The logo for SHENTEK, with the word in a bold, sans-serif font. The 'S' and 'H' are blue, while the 'E', 'N', 'T', 'E', and 'K' are green.

Huzhou Shenke Biotechnology Co., Ltd.

[www.shentekbio.com](http://www.shentekbio.com)

Address: 8th Floor, 6B Building, No.1366 Hongfeng Road, Huzhou 313000, Zhejiang Province, China

E-mail: [info@shentekbio.com](mailto:info@shentekbio.com)

Phone: +1 (908) 822-3199 / (+86) 400-878-2189